



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LIFEPAK 15 Bags (set) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Set of bags with attachments designed to be mounted on the LIFEPAK 15 defibrillator,
- Bags equipped with zippers,
- Robust material and resistant to numerous abrasions,
- Includes:
 - Main bag with storage compartments on both sides along with sturdy plastic attachments,
 - Shoulder strap for transporting all equipment,
 - Upper elongated bag for storing hard spoons operating on a zipper system,
- Dimensions:
 - Main double-sided bag: 22 x 64 x 27 cm,
 - Side pockets: 11 x 22 x 26 cm,
 - Rear pocket: 31 x 8 x 17 cm,
 - Hard spoon bag: 34 x 11 x 8.5 cm,
- Weight: 3.5 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

