



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED ŻYWIEC SU-04 Operating O/R table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Emergency power supply: through 2 new built-in 12 V / 12 Ah batteries,
- Frequency: 50/60 Hz,
- Upholstery with medical certificate in black (upholstery can be exchanged for a new one at additional cost),
- 5 - segmented,
- FAMED ŻYWIEC SU-04 operating table designed for procedures and surgeries in general surgery, vascular surgery, cardiac surgery, neurosurgery, urology, gynecology, proctology, laryngology, ophthalmology, endoscopy, laparoscopy, trauma surgery, oncology, dental surgery, plastic surgery and others,
- Electric adjustment by wired remote control:
 - Height range: 79 - 104 cm,
 - Trendelenburg position: 25°,
 - Anti-Trendelenburg position: 25°,
 - Lateral tilt in range: 20°,
 - "0" position,
- Manual adjustment:
 - Kidney uplift (kidney bench) by crank within 16 cm,
 - Infinitely adjustable backrest in the range: -40° to +70°,
 - Infinitely adjustable headrest in the range: -20° to +50°,
 - Divided side-reversible footrests with smooth adjustment in the range: -90° to +90°,
 - Angle of the split footrests: 180°,
 - Longitudinal movement of the table top mechanical in the range of 25 cm,
- Structure made of stainless steel,
- Adapted for C-arm,
- Features X-RAY translucency,
- Rails for attaching additional accessories,
- Mobile with central brake,
- The table top is removable adapted for transport by means of a dedicated cart, which makes it possible to transport the patient both to the operating room and to the postoperative room on the table top set on a transport cart (no cart included),
- It is possible to retrofit the table (for an extra charge) with additional accessories, examples shown in the last photos,
- Included:
 - Power cord,
 - Adjusting crank,
- Dimensions:
 - Bed: 200 x 50 cm,
 - Overall: 203 x 56 x 79/104 cm,
- Maximum load capacity: 200 KG,
- Weight: 265 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

