



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BIOBASE BKM-Z45B Class B laboratory autoclave (Ex-Display)

- New (after sales),
- Real photos of the product,
- Made in the PRC, (2023),
- Power supply: 380 V,
- Frequency: 50/60 Hz,
- Power: 5800 W,
- Features:
 - LCD touch screen displaying information on temperature, pressure, time, operating status, errors and alarms,
 - Automatic overheat protection, multiple evaporator control and protection, door safety lock, double overpressure protection,
 - Mobile temperature probe can directly measure the internal temperature of the liquid to ensure the sterilization effect of the liquid,
 - The sterilizer has an automatic door design, which is quick and convenient to operate, effectively protect and prevent burns,
 - High-efficiency vacuum extraction system, low noise, extreme vacuum over -90KPa, high extraction speed, good dryness, more complete sterilization,
 - Heating, sterilization, extraction and drying are completely automatic,
 - The chamber is heated by a film heater, saving energy and time,
 - High efficiency 0.22 μ m air filter to prevent recontamination of sterilized materials, with filtration efficiency of more than 99.5%,
 - The water filter in the tank and chamber can be quickly removed and easily cleaned,
 - Built-in stainless steel water tank, internal steam and water cycle,
 - Standard built-in printer,
 - Standard USB port,
- Programs:
 - Exposed Procedures,
 - Packaging procedures,
 - Rubber procedures,
 - Dressing procedures,
 - Liquid procedures ,
 - Liquid personalization,
 - Substrate sterilization procedures,
 - Substrate dissolution procedures,
 - Permanent personalization,
 - B&D tests,
 - Vacuum tests,
- Technical specifications:
 - Capacity: 45 l,
 - Chamber dimensions: ϕ 31.6 x 62.1 cm,
 - Design pressure: -0.1 / 0.3 MPa,
 - Design temperature: 144°C,
 - Temperature range: 105°C to 138°C,
- Included:
 - 3 x tray size: 24 x 46 cm,
 - Tray holder,
 - Spare door seal,
 - Printer paper,
- Dimensions: 64 x 83 x 55 cm,
- Weight: 140 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

