



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KONKRET ONYX OYK-02-1 Operating O/R table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Poland, (2011),
- Power supply: 230 V,
- Emergency power supply: via 2 built-in 12 V / 12 Ah batteries,
- Frequency: 50/60 Hz,
- Upholstery with medical approval in black,
- ONYX OYK-02-1 is a comfortable operating table, equipped with a vehicle system including free electrical and manual adjustment of each table segment,
- 5 - segment,
- Electromechanical drive,
- Electrical adjustment with a manual remote control:
 - Height in range: 68 - 115 cm,
 - Trendelenburg and anti-Trendelenburg position in range: 30°,
 - Lateral tilt: 22° / 0° / 22°,
 - Back rest: +73° / 0° / -41°,
 - Button for '0' position,
- Electrical adjustment by means of the built-in panel in the table column:
 - Height in range: 68 - 115 cm,
 - Trendelenburg and anti-Trendelenburg position in range: 30°,
 - Lateral tilt: 22° / 0° / 22°,
 - Back rest: +73° / 0° / -41°,
 - 0' position button,
- Manual adjustment by foot lever with position selection:
 - Height in range: 68 - 115 cm,
 - Trendelenburg and anti-Trendelenburg position in range: 30°,
 - Lateral tilt: 22° / 0° / 22°,
 - Back rest: +73° / 0° / -41°,
- Manual adjustment:
 - Headrest: +30° / 0° / -36°,
 - Foot segment: -90° / 0° / +28° (free side-to-side adjustment),
 - Longitudinal extension of the table in the range: 30 cm,
 - 2 - segment lumbar section via crank:
 - Push-out in range: 17 cm,
 - Lower section inclination: 42° / 0°,
 - Upper segment inclination: 22° / 0°,
- Included:
 - Wired remote control,
 - Adjustment crank,
 - Power cord,
- Dimensions:
 - Headrest: 49 x 27 cm,
 - 2-segment backrest: 49 x 55 cm,
 - Seat: 49 x 59 cm,
 - Leg segment: 23 x 64 cm,
 - Overall: 58 x 210 x 68/115 cm,
- Weight: 180 KG,
- Has an up-to-date inspection and is ready for use,
- Technical Passport (Service Report) issued, valid for a period of 12 months,
- Warranty:
 - 6-months for the domestic market (Poland),
 - 3 months for the international market,
 - Possibility to extend the warranty to 12 months for the domestic and international market for a surcharge,
- Financing options (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU

