



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ALVO MEDICAL PRELUDIUM 4-01 Operating O/R table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland, (2013),
- Upholstery with medical approval in black color,
- ALVO Preludium operating table allows safe positioning of the patient for most surgical procedures,
- With a wide range of available accessories, the table can be configured according to the requirements of different surgical disciplines,
- An excellent choice for the following disciplines: colorectal surgery, general surgery, gynecology, ophthalmology, otorhinolaryngology, urology, plastic surgery, outpatient clinics, clinics, ECPW endoscopy labs,
- 5-segment,
- Table top sections are X-RAY permeable,
- Manual adjustment:
  - Headrest via gas spring in the range:  $\pm 45^\circ$ ,
  - Back segment by gas springs in the range:  $+ 85^\circ \div - 45^\circ$ ,
  - Foot segment adjustable via gas springs range:  $+ 20^\circ \div -90$ , split and side-reversible max.  $180^\circ$ ,
  - Lateral tilt via crank in range:  $\pm 25^\circ$ ,
  - Trendelenburg and anti-Trendelenburg positions in the range:  $\pm 30^\circ$ ,
  - Table height by foot lever in range: 75-100 cm,
- Table made of stainless, acid-resistant steel,
- Equipped with:
  - Rails for attaching additional accessories,
  - Mobile base on casters,
  - Manual central brake,
  - Crank for tilt adjustment,
- It is possible to retrofit the table (for an extra charge) with additional accessories, examples shown in the last photos,
- Dimensions:
  - Overall: 55 x 210 x 75/100 cm,
  - Base: 57 x 98 cm,
- Maximum lifting capacity: 200 KG,
- Weight: 220 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

