



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **TOMEY RC-800 Auto Kerato-Refractometer (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Japan, (2015),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 100 VA,
- **Selectable menu in language: English,**
- **The RC-800 autokeratorefractometer** performs high-quality measurement combined with speed of examination and simplicity of operation, making work with **the RC-800** fast and professional,
- The device can be operated and adjusted using a joystick and touch screen, all in an extremely short time,
- All test parameters can be set on the 5.7" touch screen,
- Operating **the RC-800** is very easy. Once the optical head is adjusted to the patient's eye, the measurement is performed automatically,
- **Technical specifications:**
  - **Refraction:**
    - Contact lens measurement option,
    - PD of the patient - 50 to 86 mm,
    - VD of the patient - 0 - 16 mm,
    - Cylinders (CYL) - 0 to  $\pm 10$  D,
    - Form of cylinders +/-mix,
    - Sphere (SPH) - from -25 to +22 D,
    - Axis - from 1° to 180°,
    - Minimum pupil size - 2 mm,
  - **Keratometer:**
    - Radius of curvature - 5 to 11 mm,
    - Corneal spherical power - 30.68 to 67.50 D,
    - CYL - 0 to  $\pm 10$  D,
    - Axis - from 0° to 180°,
    - Corneal diameter measurement,
  - Has a built-in thermal printer,
  - Output port: RS232C, USB,
- **Device placed on a mobile electric table:**
  - Height adjustment in range: 73 - 98 cm,
  - Built-in 2 sockets for powering additional equipment (230 V),
  - Mobile table on 4 casters (2 braking),
- **Includes:**
  - User manual in Polish (PDF),
  - Power cord,
- Dimensions:
  - Device: 30 x 50 x 49 cm,
  - Table: 66 x 45 x 77/98 cm,
- Weight: 17 KG,
- **Has a current inspection and is ready for use,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Warranty:**
  - **6-months for domestic market (Poland),**
  - **3-months for the international market,**
  - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

