



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **YESON E-Series 8l Black LED Autoclave (New)**

New,

Modern Yeson E-series autoclave (LED) with a capacity of 8 liters,

European class B in accordance with the EN 13060 standard,

The device is designed for use in: hospitals, medical facilities (in clinics), including medical and dental offices, as well as laboratories. The autoclave can also be used for non-medical and veterinary applications.

Features and advantages of the E-series autoclave (LED):

Stainless steel chamber,

Built-in thermal printer,

Built-in printer,

Safety valve,

LED display,

Bacteriological filter,

Easy operation,

4 ready to use programs and 2 test programs,

Sterilization of packaged and unwrapped items,

Option to adjust sterilization time,

Test programs - Bowie & Dick/Helix Test, Vacuum test (leak test),

Option to maintain temperature until the next cycle,

Technical data:

Medical device: class IIb,

Chamber capacity: 8 liters,

Internal tanks for distilled and used water,

Loading method: basket with 3 trays,

Power supply: 230V  $\pm$  10% , 50/60 Hz,

Rated power: 1950 W,

Maximum working pressure: 0.21 ~ 0.24 MPa / 2.1 bar ~ 2.4 bar,

Maximum operating temperature: 134°C ~ 137°C,

Dimensions with door closed: 628.22 mm x 436 mm x 362.3 mm,

Chamber dimensions: 180 x 320 mm,

Tray dimensions: top: 13 x 27 cm | middle: 15 x 27 cm | bottom: 13 x 27 cm,

Certificate of compliance with Directive 93/42/EEC (MDD),

Net weight: approximately 43 kg

Has a current Technical Passport issued,

Warranty 24 months,

