



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET ALPHASTAR 1132 Operating O/R table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in 12 V / 17 Ah batteries,
- Frequency: 50/60 Hz,
- Power: 400 VA,
- MAQUET ALPHASTAR is a mobile universal operating table with electrohydraulic drive. Emergency control panel on the table column,
- New upholstery (memory visco) with medical approval in black,
- Adjustable in many planes by wireless remote control or through the panel on the column:
  - ON/OFF,
  - LOCK/UNLOCK - lock and unlock the table,
  - NORMAL/REV - orientation of the patient on the operating table,
  - TREND/REV - change of the patient's lying direction in relation to the table column: Trendelenburg and anti-Trendelenburg positions within  $\pm 30^\circ$ ,
  - The table is equipped with an auxiliary drive controlled by a hand-held remote control possible forward and backward travel,
  - UP/DOWN - operating surface height (up/down),
  - TILT L/TILT R - lateral tilt within  $20^\circ$ ,
  - BACK UP/BACK DN - back segment position adjustment,
  - LEG UP/LEG DN - adjustment of leg segment: up/down,
  - FLEX/REFLEX - bending the top in half in both directions: roof position / V position (Beach chair),
  - LEVEL - resetting: return to horizontal recliner position,
- Included:
  - Power cord,
  - Wired remote control,
- Dimensions:
  - Headrest: 50 x 26 x 7 cm,
  - Back rest: 55 x 54 x 7 cm,
  - Seat: 55 x 42 x 7 cm,
  - Main leg rest segment: 55 x 52 x 7 cm,
  - Lower leg rest segment: 54 x 26 x 7 cm,
- Maximum load capacity: 450 KG,
- Weight: 325 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU

