



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET ALPHAMAXX operating table 1133 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Emergency power supply: via 2 new built-in 12 V / 10 Ah batteries,
Frequency: 50/60 Hz,
Power: 400 VA,
MAQUET ALPHAMAXX
is an operating table developed in close cooperation with physiotherapists, surgeons and nursing staff. It simplifies daily work in the operating room, supports patient safety and comfort. It can be easily modified according to the needs of any surgical discipline in just a few quick steps,
New upholstery (memory visco) with medical approval in black,
Adjustable in multiple planes with a wired remote control or via a panel on the column:
Height adjustment in the range: 60 - 105 cm,
Longitudinal movement in range: 22 cm,
Zeroing - return to the horizontal position of the tabletop,
Flex and Reflex - bending the tabletop in half in both directions,
Side tilt right and left in range: 20°,
Vertical tilt - Trendelenburg and Anti-Trendelenburg position: 30°,
Backrest section adjustment: -40 / +80°,
Possibility to move the table forward / backward with electric drive,
Electric brake,
Manual adjustment:
Headrest: -45 / +90°,
Leg section - up/down -90°/10°,
Leg section - lateral retraction: 75°,
Specifications:
5 segments,
Electrohydraulic drive,
Universal table top,
X - RAY translucency,
Base housing made of durable plastic,
Equipped with rails for attaching additional accessories,
Table adapted for C-arm,
Mobile, on double casters,
Includes:
User manual in English (PDF),
Power cord,
Remote control,
Dimensions:
Headrest: 50 x 30 x 6 cm,
Upper back support segment: 52 x 24 x 6 cm,
Middle seat-back segment: 52 x 94 x 6 cm,
Foot segment of upper backrest: 24 x 12 x 6 cm,
Leg rest segment: 24 x 64 x 6 cm,
Overall length: 232 cm,
Overall width: 57 cm,
Base dimensions: 52 x 105 x 22 cm,
Maximum load capacity: 450 kg (patient + accessories),
Weight: 340 KG,
Cat 04,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

