



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED ŻYWIEC SU-10.0 Operating O/R table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Emergency power supply: through 2 new built-in 12 V / 12 Ah batteries,
- Frequency: 50/60 Hz,
- Power: 250 W,
- Upholstery (memory visco) with medical certificate in black,
- 5 - segmented,
- Made of stainless acid-resistant steel,
- Electro-hydraulic drive system ensures safe and trouble-free operation when moving the bed,
- The design of the table allows static load drainage through anti-static wheels,
- Small base facilitates surgeons' access to the patient,
- Lack of transverse metal beams and asymmetrically placed column allows intraoperative use of C-arm without the need for longitudinal movement of the tabletop,
- Change table orientation with a single button on the wired remote control,
- Display indicates battery level, table orientation and angle of each segment as it changes position,
- Infinitely variable electric adjustment of the table using the wired remote control:
 - Change of the patient's lying direction in relation to the table column: Trendelenburg position within 15° and anti-Trendelenburg position within 20°,
 - Adjusting the height of the operating surface (up / down),
 - Change of lateral tilt in the range of 18°,
 - Height adjustment in the range: 80 - 118 cm,
 - Segment:
 - Back in the range of (- 30)° / 0° / 58°,
 - Legs in the range of (- 85)° / 0° / 19°,
 - Zero function: return to horizontal recliner position,
- Possibility of manual adjustment:
 - Angle of the headrest in the basic position in the range of (- 27)° / 0° / 49°,
 - Adjustment of both parts of the split leg segment laterally within 90°,
- The table has:
 - A wide range of adjustment possibilities thanks to the option to swap all table segments in place except the fixed back segment,
 - X - RAY translucency,
 - Multifunctional wired remote control with LCD display,
 - Rails for attaching additional accessories,
 - The table is equipped with a dedicated FAMED ŻYWIEC WB-01.3 mobile cart for free movement of the table,
- Includes:
 - Power cord,
 - Wired remote control,
 - User manual in Polish,
- Dimensions:
 - Lying surface: 205 x 49 x 5.5 cm,
 - Overall: 205 x 56 x 60/84 cm,
- Maximum lifting capacity: 180 KG,
- Weight: 260 KG,
- Possibility to retrofit the table with the accessories shown in the last 2 photos (for an extra charge),
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

