



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **KARL STORZ 29003BH Endoscope cart (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- KARL STORZ 29003BH is a stable and sturdy endoscope cart with free movement,
- It has:
  - Built-in transformer,
  - Top and middle table top with dimensions: 43 x 48 cm,
  - Lower top with dimensions: 50 x 48 cm,
  - Pull-out drawer with dimensions: 39 x 38 cm,
  - Tripod with free height adjustment,
  - 3 built-in power cords for additional accessories,
  - Rubber bumpers,
- Mobile cart on 4 wheels (2 braking),
- Includes:
  - Power cord,
  - 3 power cords additional accessories,
  - 2x key for pull-out drawer,
  - Pull-out tripod,
  - Inbus key,
- Overall dimensions: 54 x 66 x 138 cm,
- Maximum load capacity:
  - Top shelf and middle shelf: 50 KG,
  - Lower shelf: 225 KG,
  - Total weight: 38 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

