



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



### **OLYMPUS CLV-180 light source (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Japanese production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 500 VA,  
Technical specifications:  
Adapted for flexible endoscopes,  
Main lighting: xenon bulb MD 631 / 300 W,  
Emergency lighting: 12 V / 35 W halogen bulb,  
Foiled control panel,  
Built-in air pump with electronic power control,  
Included:  
User manual in English (PDF),  
Power cord,  
Dimensions: 54 x 39 x 16 cm,  
Weight: 15 KG,  
Has an up-to-date inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

