



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Mobile surgical lamp COLD LIGHT YD01-5 (LED) (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
P.R.C. production, (2021),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 20 VA,  
YD01-5 lamp  
is designed for medium and small surgical operations. Easy to operate with free canopy height adjustment from 106 d0 174 cm,  
Free intensity adjustment with built-in knob,  
Adjustable with a knob to tilt the bulbs in the canopy,  
Freely adjustable rotation of the canopy and its tilt,  
Mobile on 4 casters (2 with brake),  
Technical specifications:  
Maximum light intensity: 80,000 lux,  
Color temperature 5000±500K,  
Depth of light intensity: 500 mm,  
Total illumination 275 W/m<sup>2</sup>  
,  
Included:  
Power cord,  
2 pcs. - autoclavable canopy holders,  
Dimensions:  
Canopy diameter: ø49 cm,  
Total: 48 x 100 x 195 cm,,  
Total weight: 50 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,  
The warranty period does not include bulbs and other light sources installed in the device,  
Possibility of financing (Poland only): Installments, Leasing, Loan,

