



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Infraasonic liposuction N.I.L. EVA SP2 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Belgian production, (2021),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
N.I.L. EVA SP2 system  
created for those users who have a system for INFILTRATION (infiltration pump) -.  
no pump included (can be retrofitted at additional cost)

EVA SP2 system:

Lipomatic 4 head

(can be retrofitted with the new EVAMATIC 5 head at an additional cost),

Control system,

Aspiration (suction) system,

Cannula set:

Mercedes 4 x 300,

3G Accelerator 3 x 250,

3G Accelerator 3 x 150,

1x Viterbo 4x 250,

Infiltration 2.5 x 200,

Infiltration 2.5 x 300,

Accessories:

5x pressure hose,

5x Medivac Flex 3L suction cartridge,

Suction cartridge,

2x jar,

Adjustment key,

Cleaning key,

Lipomatic ® oil (125ml),

Electric power cord,

2x Pedal,

Optional equipment (at an extra charge):

Compressor DK50 2V Silent S DRY Monsun - with soundproofing enclosure and dehumidifier,

Compressor DK50 2V DRY Monsun - without housing with dehumidifier,

MEDELA Fat Transfer Kit (Fat Collector / Fat Transfer),

EUROMI NLF1 Fat Transfer Kit (Fat Collector / Fat Transfer),

EUROMI NLF2 Fat Transfer Kit (Fat Collector / Fat Transfer),

Consumables,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

