



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE MD15T cardiac monitor with SOLAR 8000M modular monitoring system (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Monitoring system
GE SOLAR 8000M
with monitor
GE MD15T
serves as the center of the monitoring system, user controls include a keypad and control knob for easy navigation,
Modular patient monitoring system enhances clinical decision-making with a continuous flow of high-quality information. It provides robust clinical standards, including basic vital signs and advanced cardiac care,
Transmission of patient information is via the
TRAM-RAC2,
Technical specifications:
TRAM-RAC 2 module:
ECG,
TEMP/CO,
2x BP,
BP/BP,
NBP,
DISPLAY
Defib. Sync,
SPO2,
Wired modules:
ECG - Factory new 3-lead ECG cable,
SOP
2
- Factory new sensor,
NIBP - Cuff + air line,
Includes:
Power cord,
ECG - Factory new 3 lead ECG cable,
SOP
2
- Factory new sensor,
NIBP - Cuff + air line,
Connection cable set,
Power supply,
Dimensions:
Overall: 35 x 40 x 40 cm,
Total weight: 20 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

