



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OPTOPOL VC-102 colposcope (video track) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production (2015),
Power supply: 230 V,
Frequency: 50/60 Hz,
Ergonomic head design
OPTOPOL VC-102
together with the straight-arm and mobile base guarantee easy and convenient operation,
Integral glare reduction system completely eliminates reflections arising on wet surfaces,
Integrated monitor
NEOVO LCD
displays the examination "live",
Technical specifications:
OPTOPOL VC-102 head:
Freely adjustable magnification: 4x to 50x,
Working distance range: 10 - 30 cm,
Field of view range: 6 - 75 mm,
Dual-zone ring illuminator composed of ultra-bright LEDs,
Image retention function,
Autofocus function or manual adjustment,
RED FREE function - electronic infrared or negative filter,
SCALE function - magnification marker or object size,
Glare reduction,
Possibility to connect ultrasound,
Green negative filter,
Working range of the arm: 60 cm,
NEOVO X-15AV monitor:
Screen diagonal: 15" inches,
Resolution: 1024 x 768 px,
TN matrix,
Format: 4:3,
Brightness: 400 cd/m²
,
Includes:
Power cord,
Overall dimensions: 88 x 70 x 134 cm,
Weight: 45 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
The warranty period does not include bulbs and other light sources installed in the device,
Possibility of financing (Poland only): Installments, Leasing, Loan

