



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Class B autoclave GETINGE K7+ (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Polish production, (2017)

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 2200 W,

Autoclave

GETTING K7+

ideally meets the sterilization needs of clinics from small dental offices to hand-held sterilization in surgical departments,

Extremely fast, the Class B cycle takes 25 minutes including drying, the Flash program takes only 9 minutes, and has a special program for

sterilizing dental handpieces, which takes 15 minutes for 10 packed handpieces,

Programs:

Class B cycle, instruments in packages, 134°C, process time ~25 min,

Class B cycle, packaged instruments, 121°C, process time ~39 min,

Class B cycle, textile materials, 134°C, process time ~29 min,

Instant program, 134°C, process time ~9 min,

Dental special program, 134°C, process time for 10 packaged tips ~15 min,

Class B cycle, prions, 134°C, process time ~38 min,

Test programs: Bowie & Dick, leak test,

Technical specifications:

Sterilization of batches type: A and B,

Chamber capacity: 20 l,

Maximum pressure: 2.4 bar,

Clean water tank: 5 l.,

Water consumption per cycle: 0.2 - 0.4 l.,

Menu in language: Polish,

Included:

Power cord,

5 x tray with dimensions 18.5 x 42 cm,

Tray holder,

Dimensions: 57 x 52 x 44 cm,

Weight: 80 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

