



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS INTELLIVUE MP30 NEONATAL patient monitor/cardiac (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Emergency power supply: through 2 built-in batteries 10.8 V / 8.7 Ah
(one included - can be retrofitted with a second new one at additional cost)

Frequency: 50/60 Hz,

IntelliVue MP30

NEONATAL

is a multifunctional cardiac monitor designed for adults, children and neonates,

It has:

Panel screen, made of TFT LCD technology with a diagonal of 10.4" inches and working in SVGA resolution,

Touch screen operation, and navigation joystick, this allows direct access to many functions with simple commands,

Convenient handle and sturdy case make it easy to carry,

Trends of measurements can be presented graphically, combining up to three curves on a single graph, making it easy to identify changes in the patient's physiological state,

Trends of fast-variable parameters can be updated every single beat and up to four high-resolution trend segments can be displayed,

IntelliVue technology - access information from various systems in the hospital network - pharmacy, laboratory, ECG management, radiology department PACS, server-based hospital systems and others - and deliver it to the patient monitor,

Lacking a separate hard drive and fan, IntelliVue monitors are small, quiet, stable and quick to boot up,

A navigation dial makes it easy to enter information and navigate the screen,

Audible and visual alarms,

LAN,

Equipped with PHILIP Model M3001A module with wiring:

ECG /RESP - New 3 - lead cable,

SOP

2

- Factory new sensor,

NBP - New cuff with air tubing,

TEMP - Temperature connector socket,

(sensor not included, can be retrofitted with a new one at additional cost),

Menu in language: Polish,

Included:

User manual in English (PDF),

Power cord,

New 3-lead ECG cable,

Factory new NIBP cuff,

Factory new SOP sensor

2

Overall dimensions: 21 x 34 x 28 cm,

Weight: 6.5 KG,

Cat 03,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

