



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **YONKER YK-8000C patient monitor/cardiac monitor (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
P.R.C. production, (2020),  
Power supply: 230 V,  
Emergency power supply: via built-in 11.1 V / 2200 mAh battery,  
Frequency: 50/60 Hz,  
Power: 70 VA,  
Advanced multi-parameter cardiac monitor  
YONKER YK-8000C  
, which can be used during surgery, post-surgery recovery, post-anesthesia and rescue to monitor the vital signs of adults, children and newborns. YK-8000C is a handy - ergonomic, comfortable and intuitive to use device,  
Cardiomonitor designed for adults, children and newborns. Age group selection in the settings with automatic adjustment of alarm limits,  
Technical specifications:  
Screen: 12.1" inch LCL TFT,  
Built-in measurement modules:  
ECG,  
SPO2,  
TEMP1 | TEMP2,  
NIBP,  
The cardiac monitor works with a central monitoring station,  
Adjustable alarm limits,  
Memory of measured values (tabular and graphical trends),  
Included:  
Power cord,  
ECG 5-lead cable,  
SPO2 sensor,  
NIBP cuff with air tubing,  
Dimensions: 30.5 x 16.2 x 29 cm,  
Weight: 3.3 KG,  
Cat 01,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

