



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Patient monitor/cardio monitor YONKER YK-8000C (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
Produced by P.R.C., (2019),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 70 VA,  
Advanced multi-parameter cardio monitor YONKER YK-8000C, which can be used during surgery, post-surgery recovery, after anesthesia and rescue to monitor vital signs of adults, children and newborns. YK-8000C is a handy - ergonomic, comfortable and intuitive to use device,  
Cardio monitor intended for adults, children and newborns. Age group selection in settings with automatic adjustment of alarm limits,  
Technical specifications:  
Screen: 12.1" inch LCL TFT,  
Built-in measurement modules:  
ECG,  
SPO2,  
TEMP1 | TEMP2,  
NIBP,  
The cardiomonitor cooperates with the central monitoring station,  
Alarm limit adjustment,  
Memory of measured values (tabular and graphic trends),  
Included:  
Power cord,  
5-lead ECG cable,  
SPO2 sensor,  
NIBP cuff with air tube,  
Dimensions: 30.5 x 16.2 x 29 cm,  
Weight: 3.3 KG,  
Cat 02,  
It has a current inspection and is ready to work,  
Technical Passport (Service Report) issued valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,  
Possibility of financing (Poland only): Installments, Leasing, Loan,

