



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



ASPEL ASCARD A4 MODEL D ECG Apparatus (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production,
Power supply: 230 V,
Emergency power supply: via 1 new built-in 12 V / 1.3 Ah battery,
Frequency: 50/60 Hz,
Power: 20 VA,
Specifications:
ECG signal: 12 standard leads,
Printout: 3, 6, 12 (after automatic analysis) - real-time channel (paper width: 112 mm),
Paper: thermo-sensitive, dust-free with a width of 112 mm,
Recording speed: 5 / 10 / 25 / 50 mm/s,
Test sensitivity: 2.5 / 5 / 10 / 20 mm/mV,
Digital filters: 50 Hz and 35 Hz; Anti-dry filter,
Built-in LCD screen: 240 x 128 characters,
Foiled, touch and readable control panel,
Functional and alphanumeric keypad,
Signaling of bad connection of individual electrodes,
Output: RS-232,
Selectable menus in the following languages: PL, ENG, GER, FR, IT, ESP, HR, TR, CS, RUM, SK,
Included:
User manual in Polish (PDF),
10-lead ECG cable,
EKK clip-on limb electrodes - 4 pcs,
EPP suction cup electrodes - 6 pcs,
Power cord,
Dimensions: 31 x 23 x 7 cm,
Weight: 3 KG,
Cat 03,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,



