



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **NORMATEC 3 HYPERICE leg lymphatic drainage system (Ex-Display)**

New (post-show),  
Real photos of the product,  
Production , (2023),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Massage and lymphatic drainage system  
NORMATEC 3  
Uses state-of-the-art compression technology to quickly regenerate tired and achy legs,  
Hyperice  
offers 7 levels of compression and patented Pulse technology for precise massage in 5 overlapping zones. In addition,  
NORMATEC 3  
is ZoneBoost™ technology that uses biomimicry to recreate natural circulation in the legs,  
Design of the drainage system  
NORAMTEC 3  
features an intuitive simple interface that allows us to easily manage the pressure therapy session before and during treatment. Bluetooth connectivity also allows the Normatec device to be managed via the Hyperice app connected to your phone,  
Features:  
Accelerate muscle recovery,  
Support in improving lymphatic drainage,  
Improved circulation,  
Relief of leg pain,  
Firms and improves skin elasticity,  
Reduction of inflammation,  
Increase muscle fitness and performance,  
HYPERICE application:  
Normatec 3 can connect to your phone and collect data via the app to help customize your massage,  
Ability to safely and easily unlock the next levels and customize the Normatec 3's action to suit your own needs,  
Advantages:  
Comfortable adjustment of time and pressure intensity settings,  
Ability to activate the ZONE BOOST function (zone amplification),  
Ability to save and easily restore preferred settings,  
Monitoring and analysis of individual statistics to optimize recovery efficiency,  
Specifications:  
7 compression levels and 5 air chambers (customizable),  
Intuitive control unit,  
Lightweight control unit (1.45 KG),  
Zone Boost™ function,  
Pressure range: 30-110 mmHg,  
Connectivity to Hyperice app,  
Operating time: up to 3 hours on a single charge,  
Durable materials and top quality workmanship,  
Legs suitable for people up to 161 cm tall,  
Included:  
NORMATEC control unit 3,  
Short version legs,  
Pressure cables,  
Power supply,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

