



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **POLDERMA ML BEAUTY MACHINE diode laser 808-A2 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Polish production, (2022).  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Diode laser is one of the most modern and safest methods of permanent hair removal. The principle of diode laser is based on the phenomenon of selective thermolysis. The laser light is selectively absorbed by the melanin contained in the hair, then turns into thermal energy destroying the hair structures responsible for its growth. The destruction of the hair follicle occurs, while at the same time the tissue surrounding it remains intact,  
The device uses the most effective wavelength.  
:  
755nm - Wave well and selectively absorbed by melanin. It works slightly shallower than 808nm so it is less painful,  
808nm - The wave is also well absorbed by melanin,  
1064nm - A wave that acts deeper than 755nm,  
The device uses an innovative 12-bar emitter,  
The laser provides high efficacy against dark and light hair regardless of phototype,  
Advantages:  
Negligible soreness,  
Permanent depilation results,  
Effective and reliable cooling system (-10 to +3 Degrees Celsius),  
Large, easy-to-read LCD display,  
Technical specifications:  
Type: diode laser,  
Wavelength: 808 / 810nm (12Bar),  
Power consumption: 2000W,  
Output power: 800W,  
Screen: 8.4" color touch screen,  
Pulse length: 10-200ms with smooth adjustment,  
Energy: 1-150 J/cm  
2  
adjustment,  
Frequency: 1-10Hz,  
Cooling: semiconductor, water, air (-10 to +3 degrees Celsius),  
Selectable menu in: English, Polish,  
Included:  
User manual in Polish,  
Treatment head,  
Dedicated box for transporting the device,  
2x box key,  
Head holder,  
Power cord,  
Dimensions:  
Device: 59 x 40 x 36 cm,  
Weight: 40 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

