



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



PENTAX FI-16RBS intubation bronchofiberscope (Ex-Display)

New (post-show),
Real photos of the product,
Japanese manufacture, (2017),
Pentax FI-16RBS
is a wireless intubation bronchofiberscope ideal for emergency use,
It provides maximum comfort,
Tilting tip makes it easy to orient and position endotracheal tubes,
Ergonomic, lightweight body provides easy control and productivity,
Equipped with CSV (Closed Suction Valve) mechanism for accurate suction,
Bright and clear high-resolution image for excellent observation,
Can be connected to a conventional light source,
Specifications:
Optical system:
Field of view: 95°,
Depth of field: 3 - 50 mm,
Endoscopic probe:
Outer diameter of endoscopic probe tip: 5.1 mm,
Outer diameter of endoscopic probe tube: 5.2 mm,
Working length: 600 mm,
Biopsy channel:
Inside diameter of the channel: 2.55 mm,
Tip bendable in range:
160° (up),
130° (down),
Included:
Led Test adapter
used to check battery life,
Silicone oil
PENTAX OF-Z11
,
One-sided brush
PENTAX CS-C3S
,
Suction control valve
PENTAX OF-B53
,
Cap
PENTAX BS-B1
,
Cylinder
PENTAX OF-B103
,
Reusable inlet valve
PENTAX OF-B190
,
Bite blocker
PETNAX OF-Z5
,
Dedicated PENTAX case for transporting instrumentation,
2 pcs key for the case included,
Dimensions:
91.5 x 32.5 x 11 cm,
Weight: 4 KG,
Cat 02,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
The warranty period does not include bulbs and other light sources installed in the device,
Possibility of financing (Poland only): Installments, Leasing, Loan,

