



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rigid endoscope (laparoscope) RICHARD WOLF PANOVIEW 8934.432 (Ex-Display)

New (post-show),
Realistic pictures of the product,
German manufacture,
Technical specifications:
Diameter: 10 mm,
Viewing angle: 30°,
Autoclavable,
Built-in fiber optic light transmission,
Includes:
Dedicated container for transporting the endoscope,
Dimensions:
Working length: 30.5 cm,
Overall length: 37.5 cm,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

