



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS KV-5 endoscopic suction machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 130 VA,
Simple and intuitive operation,
Suitable for all endoscopic procedures,
Variable suction with vacuum pressure gauge: 0 to -1 bar,
Maximum suction pressure: -85 kPa,
Quiet pump operation,
Includes:
Power cord,
User manual in English (PDF),
Seeres jar with 1000 ml disposable cartridge,
Hydrophobic microbial filter,
Dimensions: 27 x 39 x 22 cm,
Weight: 13 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

