



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ECG Device Electrocardiograph WELCH ALLYN CP 200 (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Manufactured in the U.S.A.,

Power supply: 230 V,

Backup power supply: via 4 new built-in 6.4 V / 6 Ah batteries,

Frequency: 50/60 Hz,

The Welch Allyn CP 200 is a 12-lead resting ECG system that can also be used for spirometry tests. The CP 200 is equipped with an alphanumeric keyboard and an LCD display. The display can show a quick color preview of the ECG reading with the ability to edit settings. The Welch Allyn system provides full-page (8.25 x 11) ECG printouts and can store up to 50 ECG measurements and 50 spirometry measurements in its internal memory. The ECG system offers the ability to transfer patient data via SD card, UBS system or your facility's EMR system, High-resolution thermal printer generates full-size standard and custom reports,

Full alphanumeric keyboard and LCD display for quick and easy entry of patient information,

Advanced filters for optimal ECG quality ensure accurate readings,

Configurable report formats allow you to print information the way you want,

Battery operation allows easy transport of the device to the patient,

Instant-on feature allows you to quickly turn on the device and immediately perform an ECG test,

Technical specifications:

12-lead resting ECG,

LCD screen for easy viewing,

Performs spirometry tests,

Stores up to 50 ECG and 50 spirometry readings,

Full-page prints (8.25 x 11),

Configurable report formats,

Cart Welch Allyn:

Mobile on 4 wheels (2 braked),

2 side baskets for accessories,

Cable stand,

Possibility to select the menu in the following languages: EN, DE, FR, ES, IT, PT, SE, NR, NL, FI, PL, RU, GR, HR, ML,

Keyboard on the control panel in English,

Included:

User manual in English (PDF),

Power cord,

Mobile transport trolley,

10-lead ECG cable,

Dimensions:

Welch Allyn CP 200: 41 x 40 x 16 cm,

General: 55 x 60 x 102 cm,

Weight: 23 kg,

Cat 02,

It has a current inspection and is ready to work,

Issued Technical Passport (Service Report) valid for the period 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,

Financing options (Poland only): Installments, Leasing, Loan,

