



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **NOUVAG AG DISPENSER DP 30 tumescent anesthesia pump (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Swiss production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 40 VA,  
Dispenser DP 30  
is a universal pump used for tumescent anesthesia in endovascular, liposuction and other procedures,  
Intuitive to use,  
Equipped with a metal stand for hanging a bottle of tumescent fluid,  
Infinitely precise flow rate control with potentiometer from 0 - 100% ,  
Foot switch allows you to turn the flow on and off and control the discharge speed without interrupting the operation,  
Specifications:  
Flow rate: 1 - 15 L/h,  
Maximum pressure: 2 bar,  
Includes:  
38 cm high tripod for hanging the tumescent fluid bottle,  
NOUVAG Model VARIO foot attachment,  
Multi-language instruction manual (PDF),  
Power cord,  
Dimensions: 25.5 x 23 x 48 cm,  
Weight: 3.5 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

