



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS PAGE WRITER 200 electrocardiograph ECG machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V,
Emergency power supply: through 1 new built-in 6 V / 7.2 Ah battery,
Frequency: 50/60 Hz,
Electrocardiograph
PHILIPS PAGE WRITER 200
is offers 12-lead cardiographic ECG readings,
The system has a built-in printer that provides printing of full ECG readings,
The device features an alphanumeric keypad and a built-in 5" inch display,
Features:
ECG signal: 12 standard leads,
Printout: 3 / 6 / 12 - channel real-time and continuous,
Recording speed: 5 / 10 / 25 / 50 mm/s,
ACI-TIPI system helps predict the likelihood of acute cardiac ischemia,
TPI helps identify patients with acute myocardial infarction (AMI) who may benefit from thrombolytic treatment,
Standard analysis programs for adults and children,
Full-page reports in English (A format) or metric (A4 format),
Selectable menus in ENG, FR, GER, GER, ESP, IT, SWE, FI, NOR, ,
Included:
User manual in English (PDF),
Power cord,
10-lead ECG cable,
Dimensions: 39 x 44 x 11 cm,
Weight: 8.5 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

