



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DEXTA ophthalmic chair MK80DX-610-SV (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 750 VA,
New medical approved upholstery in green color,
Ophthalmic chair MK80DX-610
is a multifunctional chair ideally suited for ophthalmology, ENT surgeries,
Through its built-in features, it will also be ideal for aesthetic medicine and also for various types of procedures performed,
Technical specifications:
Manual adjustment:
Free adjustment of the headrest with knobs in extension and tilt,
Possibility to rotate the chair on both sides in the range of 45°,
Electric adjustment by foot control:
Backrest adjustment within: -6° / 0° / 72° (automatic synchronized adjustment with leg rest),
Height adjustment within: 70 - 95 cm,
Adjustment of seat tilt in Trendelenburg position range: 16°,
4 programmable chair functions,
The chair has built-in English tup rails as standard
(possibility to retrofit the chair with EUROSHINS at an additional cost),
Included:
Power cord,
Wired foot controller,
Dimensions:
Overall: 77 x 110 x 128/153cm,
Headrest: 22 x 19cm,
Back rest: 45 x 46 cm,
Seat: 52 x 50 cm,
Leg rest: 50 x 45 cm,
Total length of reclining position: 168 cm,
Weight: 110 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

