



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Scrub Sink Triple Uzumcu 80800 (New)**

- New
- Made in Turkey
- Made of 304 quality stainless steel designed to be used by the surgery team for hand and nail cleaning before entering the sterile zone
- Photocell display screen that shows the temperature and duration of water alternately from each faucet
- Specially designed with round concave corners, easy to clean and antibacterial material considering the areas of use & user's requests and to be compatible with international hygienic norms
- Photocell and knee controlled water flow and knee controlled soap flow at each station
- LCD screen displays water temperature and washing time counter
- Thermostatic mixing valve for cold and hot water
- Dimensions: 2200 x 675 x 1250 mm,
- Technical Passport
- 24 months warranty

