



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ultrasonic Cleaning Unit Uzumcu PM-55 (New)

- New
- Made in Turkey
- Ultrasonic power : 1500 W / 750 eIf
- Frequency : 28 kHz
- Heating power : 1500 W
- Sink capacity : 40 Lt
- 0 - 90°C analogue temperature adjustment
- 0 - 30 minute analogue timer
- Made of 304 quality stainless steel to be used in the decontamination areas of sterilization units of hospitals in line with user expectations
- One or two washing sinks depending on user requirements
- Ultrasonic washing section
- 3 cabinets with double doors each
- Faucets with hot & cold water connections for each sink
- Air/water gun with 8 different features
- Dimensions:
 - standard sink: 500 x 400 x 250 mm,
 - ultrasonic sink: 500 x 400 mm,
 - Total: 2700 x 600 x 900 mm,
- Technical Passport
- 24 months warranty

