



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## Ultrasonic Cleaning Unit Uzumcu PM-55 (New)

- New
- Made in Turkey
- Ultrasonic power : 1500 W / 750 eIf
- Frequency : 28 kHz
- Heating power : 1500 W
- Sink capacity : 40 Lt
- 0 - 90°C analogue temperature adjustment
- 0 - 30 minute analogue timer
- Made of 304 quality stainless steel to be used in the decontamination areas of sterilization units of hospitals in line with user expectations
- One or two washing sinks depending on user requirements
- Ultrasonic washing section
- 3 cabinets with double doors each
- Faucets with hot & cold water connections for each sink
- Air/water gun with 8 different features
- Dimensions:
  - standard sink: 500 x 400 x 250 mm,
  - ultrasonic sink: 500 x 400 mm,
  - Total: 2700 x 600 x 900 mm,
- Technical Passport
- 24 months warranty

