



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Endoscopy Workstation OLYMPUS (CLV-180, CV-180) (Reconditioned)**

- Condition: Reconditioned - Very Good,
- Made in Japan,
- Compartments:
- **Light Source OLYMPUS Model CLV-180 EVIS EXERA II:**
  - Compatible with flexible endoscopes,
  - Xenon Light,
  - Power supply: 230 V,
  - Frequency: 50 - 60 Hz,
  - Power: 500 VA,
  - Xenon bulb: MD 631 300 W,
  - Emergency light: halogen 12 V 35 W,
  - Sealed control panel, touchscreen,
  - Air pump,
  - Water container,
  - Dimensions: 38,5 x 53 x 17 cm,
- **Video System Center OLYMPUS Model CV-180 EVIS EXERA II:**
  - Power input: 230 V,
  - Frequency: 50 - 60 Hz,
  - Power: 150 VA,
  - High Quality image in HDTV
  - Compatible with: EVIS 100/130/140 and EVIS EXERA 160,
  - Output: HD / SD SDI,
  - Sealed control panel,
  - Comes with:
    - PigTail Maj cable 1430 with handle,
    - Keyboard: Olympus EVIS EXERA MAJ 1426,
    - Memory card xD with PC adapter,
  - Dimensions: 38,5 x 49 x 10 cm,
- **Gastroscope Video Endoscope OLYMPUS Model Evis Exera II GIF Q-180:**
  - Condition: Reconditioned - Very Good,
  - Made in Japan,
  - Direction of View: Forward,
  - Outer diameter: 9 mm,
  - Outer diameter insertion tube: 8,8 mm,
  - Working field diameter: 2,8 mm,
  - Working length: 1030 mm,
  - Total length: 1345 mm,
  - Depth of Field: 3-100 mm,
  - Max Angulation Up/Down: 210°/90°
  - Max Angulation Right/Left: 100°/100°,
  - Field of View: 140°,
  - Compatible with processors: CV-160/140,
  - Original suitcase: 59 x 43 x 16 cm,
  - Valves,
  - Manual in English (pdf),
- **Medical Display OLYMPUS Model OEV 191H:**
  - Made in Japan,
  - Condition: Reconditioned - Very Good,
  - LCD 19",
  - Resolution: 1280 x 1024,
  - Output: composite, S-video, RGBS, VGA,
- **Endoscopy Video Cart OLYMPUS WM-NP1:**
  - Condition: Reconditioned - Very Good,
  - Made in UK,
  - Power Input: 230 V,
  - Frequency: 50 - 60 Hz,
  - **Equipped with:**
    - 3 Articulated Arm for LCD,
    - 2 Articulated Arm with keyboard stand 40 x 21 cm,
    - Gas cylinder handle,
    - 2 shelves : 39,5 x 52 cm,
    - Countertop : 38 x 47 cm,
    - Bottom Shelf : 48 x 61 cm,
    - Electrical System,
    - Handle,
    - Endoscope handle,
  - Adjustable shelves height,
  - On casters,
  - Dimensions: 65 x 65 x 150 cm,
- Technical Passport + check-up,
- 30-days warranty

