



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Examination Table EOS E04 (New)**

New,

Made in Poland,

It has a stable steel construction coated with powder paint,

Standard features include three-section sleep surface, foldable to Pivot and seat positions, made of plywood,

Table height adjustment is carried out with an electric actuator operated by a manual remote control,

Basic option - manual remote control - free of charge,

Convertible option - selection of foot remote control (no hand control) - for extra charge,

Complete option - hand control + foot control - at extra charge,

Fully ball-bearing mounted height adjustment mechanism,

Individual parts of the tabletop are adjustable with gas springs or, as an additional option, with an electric actuator controlled with the hand-held remote control,

There is a free option to order upholstery in the colour of your choice, except for "PRESTIGE" colours which are available for an extra charge,

Table features and benefits

Gas spring by German company BANSBACH,

Electric actuator by a German company DEWERT,

Upholstery of very good quality,

5 cm thick foam,

Steel ball bearings were used in the movable parts of the table, which makes the table very reliable,

Frame elements made of the best Polish steel,

As additional options are available (for an extra charge):

Armrests,

Side shelves,

Hanger for the underlay,

Castor system,

Foot remote control (instead of manual remote control),

Control from the frame,

Euro rails (40cm),

Knee bolsters,

Injection support,

Arm support for the chair,

Electric adjustment of the table top angle to the chair position,

Electric adjustment of tabletop angle to Pivot position,

Technical parameters:

Length 200 cm,

Width 65/70/80 cm,

Height 48 - 100 cm,

Headrest angle adjustment -700 to 400,

Backrest angle adjustment 00 to 850,

Adjustment of the angle of the backrest 00 to 300,

Weight 91 kg,

Maximum load 160 kg,

It has a current Technical Passport,

Warranty on structure 24 months,

Warranty on hydraulic accessories 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPAÑOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANÇAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

