



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **RODENSTOCK CV 900 LCD panel (New)**

- New,
- Made in Germany,
- Advantages and features of the RODENSTOCK CV 900 LCD panel:
  - Infrared connection to Phoromat 2000,
  - Red-green optotypes for binaural vision testing,
  - Modern solution for ophthalmology offices,
  - Programmable distance: patient - optotype,
  - Convenience of use,
  - Wall mounting,
  - Distance sequence program,
  - Color tests,
  - Contrast sensitivity test,
  - Unlimited number of tests,
  - Patient education program,
  - IR sensor,
  - Red and green glasses included,
- Specifications:
  - Display: 22" LCD panel,
  - Distance units: feet or meters,
  - Working distance: 2 to 7 m (measuring step 0.1 m),
  - Program: 3 programs with up to 15 optotypes,
  - External input/output port: 2 x USB,
  - Communication: infrared for wireless connection,
  - Remote control: included,
  - IR remote control: 3 available IR channels,
- AVAILABLE OPTOTYPES AND TESTS:
  - Standard: Snellen, Landolt, letters, numbers, characters,
  - Color vision: Ishihara, 100 HUE,
  - Red - green: Schober, Worth, Bichrom-Balance, Fixation Disparity,
  - Poor visibility: ETDRS, Crowding Bars, AMD-Optotypes, Amsler,
  - Contrast: Contrast sensitivity, Contrast adjustment,
  - Children: HOTV, hands, images, symbols, animations,
  - Additional: Hearing/speech impaired, patient education, 3D images,
- DIMENSIONS AND ELECTRICAL REQUIREMENTS:
  - WDH dimensions: 605 x 30 x 370 mm,
  - Weight: Approx. 3.8 kg,
  - Power supply: AC 100 V to 240 V,
  - Frequency: 50/60 Hz,,
  - Power consumption: Maximum 35 W,
- Has a current Technical Passport issued,
- Warranty 12 months,

