



Examination Tray CATER (New)

- New
- Made in Poland
- Stable frame made from powder coated steel
- On castors (with breaks)
- Countertops available in 3 colors
- **Additional equipment (paid):**
 - Tempered glass shelf: 6 mm
 - Chromium-nickel shelf
 - Two handles for bowls: plastic or stainless steel,
- Technical Data:
 - Width: 78 cm,
 - Height: 80 cm,
 - Depth: 45 cm,
 - Metal countertop: 16,5 kg,
 - Glass countertop: 20 kg,
 - Max load: 15 kg,
- 12 months warranty

