



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006

kw@praiston.pl



Rehabilitation wedges KLIN (New)

New,

Polish production,

They are designed to perform exercises and rehabilitation procedures aimed at treatment or alleviation of diseases and injuries in an effective way,

Balls are made of polyurethane foam covered with high quality leather-like material,

Available sizes (length/width/height) (cm):

20x15x10

20x25x10

30x20x16

38x38x5

38x38x7

30x40x12

30x40x16

35x15x21

40x50x20

50x35x20

60x60x12

70x60x25

70x60x30

70x60x35

80x60x40

100x80x30

120x60x40

Other sizes are available on request,

12 month warranty,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

