



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



Rehabilitation wedges KLIN (New)

New,
 Polish production,
 They are designed to perform exercises and rehabilitation procedures aimed at treatment or alleviation of diseases and injuries in an effective way,

Balls are made of polyurethane foam covered with high quality leather-like material,

Available sizes (length/width/height) (cm):

- 20x15x10
- 20x25x10
- 30x20x16
- 38x38x5
- 38x38x7
- 30x40x12
- 30x40x16
- 35x15x21
- 40x50x20
- 50x35x20
- 60x60x12
- 70x60x25
- 70x60x30
- 70x60x35
- 80x60x40
- 100x80x30
- 120x60x40

Other sizes are available on request,

12 month warranty,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

