



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation wedges KLIN (New)

New,
Polish production,
They are designed to perform exercises and rehabilitation procedures aimed at treatment or alleviation of diseases and injuries in an effective way,

Balls are made of polyurethane foam covered with high quality leather-like material,

Available sizes (length/width/height) (cm):

20x15x10
20x25x10
30x20x16
38x38x5
38x38x7
30x40x12
30x40x16
35x15x21
40x50x20
50x35x20
60x60x12
70x60x25
70x60x30
70x60x35
80x60x40
100x80x30
120x60x40

Other sizes are available on request,

12 month warranty,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

