



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Aesthetic Medicine chair PROMAT MX (New)

New,
Germany production,
Power supply: 230 V,
Frequency: 50/60 Hz,
The standard color of the armchair is GRAFITOWO-WHITE,
Table equipped with wide armrests ,
Equipped with 4 electric motors,
The table's functions can be adjusted with a manual remote control (standard or at an extra charge with a foot controller),
Controlled by remote control or foot controller:
Height adjustment in range: 55 - 81 cm,
Backrest adjustment in range: 0° - 78°,
Seat adjustment in range: 0° - 33°,
Adjusting the foot of the calf section in range: 0° - 75°,
Manual adjustment:
Headrest adjustment: 180°,
Footrest adjustment heel part: 180°,
Adjustment of armrests in range: 0° - 180°,
Maximum table load capacity: 175 kg,
Retractable casters for easy manoeuvrability,
Additional equipment for an extra charge:
Handle for mounting the paper roll,
Foot controller,
Multifunctional headrest with breathing hole,
Possibility of changing the color of upholstery,
Noise level: 60 dBA,
Weight: 65 kg,
Size:
Headrest: 240 x 375 mm,
Backrest: 610 x 545 mm,
Seat: 520 x 550 mm,
Calf part footrest: 275 x 525 mm,
Footstool heel part: 195 x 445 mm,
Armrest: (2x) 470 x 75 mm,
Total length: 182 - 209,5 cm,
Total width: 54.5/74.5 cm (including armrests),
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan

