



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Chair PROMAT NG (New)

- New,
- Made in Germany,
- Color: Black & White,
- Equipped with 4 powerful actuators,
- Maximum lifting capacity up to 240 kg,
- Operated by foot switch or remote control,
- The retractable wheels are convenient for easy manoeuvring,
- **Additional Equipment (extra paid):**
 - Multifunctional headrest with breathing hole (when lying face down),
 - Paper dispenser,
 - Foot control,
 - Handle for accessories,
 - Handle for drip,
 - Handle for remote control,
 - Leg Rest assembly system,
 - Back Rest assembly system,
 - 7,2 Ah powerful actuators,
 - Choice of upholstery colours,
- Technical Passport,
- 24-month warranty

