



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Transfer Chair NITROCARE NTS X6 (New)

- New
- Made in Turkey
- Easy of use with ergonomic design, light structure
- Excellent maneuverability
- Adjustable arm supports
- Armrests with gradual opening and closing
- Central brake bar
- Adjustable headrest, backrest
- Polyurethane backrest and seating
- Dimensions;
 - Height: 1160 mm
 - Depth: 1000 mm,
 - Width: 643 mm,
 - Seat depth: 440 mm,
 - Backrest width: 500 mm,
 - Armrests: 355 x 425 x 65 mm,
- Technical Passport
- 24 months warranty

