



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Transfer Chair NITROCARE NTS X6 (New)

- New
- Made in Turkey
- Easy of use with ergonomic design, light structure
- Excellent maneuverability
- Adjustable arm supports
- Armrests with gradual opening and closing
- Central brake bar
- Adjustable headrest, backrest
- Polyurethane backrest and seating
- Dimensions;
  - Height: 1160 mm
  - Depth: 1000 mm,
  - Width: 643 mm,
  - Seat depth: 440 mm,
  - Backrest width: 500 mm,
  - Armrests: 355 x 425 x 65 mm,
- Technical Passport
- 24 months warranty

