



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Examination Table PROMAT NG (New)**

- New,
- Made in Germany,
- Color: Black & White,
- Equipped with 4 powerful actuators,
- Maximum lifting capacity up to 240 kg,
- Operated by foot switch or remote control,
- The retractable wheels are convenient for easy manoeuvring,
- **Additional Equipment (extra paid):**
  - Multifunctional headrest with breathing hole (when lying face down),
  - Paper dispenser,
  - Foot control,
  - Handle for accessories,
  - Handle for drip,
  - Handle for remote control,
  - Leg Rest assembly system,
  - Back Rest assembly system,
  - 7,2 Ah powerful actuators,
  - Choice of upholstery colours,
- Technical Passport,
- 24-month warranty

