



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **PROMAT MX treatment table (New)**

New,  
Germany production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
The standard color of the armchair is GRAFITOWO-WHITE,  
Table equipped with wide armrests ,  
Equipped with 4 electric motors,  
The table's functions can be adjusted with a manual remote control (standard or at an extra charge with a foot controller),  
Controlled by remote control or foot controller:  
Height adjustment in range: 55 - 81 cm,  
Backrest adjustment in range: 0° - 78°,  
Seat adjustment in range: 0° - 33°,  
Adjusting the foot of the calf section in range: 0° - 75°,  
Manual adjustment:  
Headrest adjustment: 180°,  
Footrest adjustment heel part: 180°,  
Adjustment of armrests in range: 0° - 180°,  
Maximum table load capacity: 175 kg,  
Retractable casters for easy manoeuvrability,  
Additional equipment for an extra charge:  
Handle for mounting the paper roll,  
Foot controller,  
Multifunctional headrest with breathing hole,  
Possibility of changing the color of upholstery,  
Noise level: 60 dBA,  
Weight: 65 kg,  
Size:  
Headrest: 240 x 375 mm,  
Backrest: 610 x 545 mm,  
Seat: 520 x 550 mm,  
Calf part footrest: 275 x 525 mm,  
Footstool heel part: 195 x 445 mm,  
Armrest: (2x) 470 x 75 mm,  
Total length: 182 - 209,5 cm,  
Total width: 54.5/74.5 cm (including armrests),  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Installments, Leasing, Loan

