



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LED Lamp LID - Bella (New)

- New,
- Made in France,
- Watts: 17 W,
- Kelvins: 4000 K,
- Lumens: 1440 Lm,
- Light source: 12 LED,
- LED lifetime: 50 000 hours minimum,
- Color Rendering Index (IRC) 95,
- Power supply cable 3 m,
- Insulation Class II,
- Available colors : white
- Dimensions :
 - Head: 160 x 133 mm,
 - Base: 69 x 69 x 115 mm,
 - Weight: 2,3 kg,
 - Length: 700 mm,
- Diameter of illuminated field :
 - Distance 30 cm: 60000 lux, illuminated field: 5 cm
 - Distance 50 cm: 45000 lux, illuminated field: 10 cm,
 - Distance 100 cm: 15000 lux, illuminated field: 20 cm,
- Double intensity switch in option (on the top of the lamp),
- The ergonomic handle allows a precise and intuitive positioning,
- Easy cleaning and disinfection thanks to the smooth surfaces and quality materials,
- LED are implanted on a printed circuit cooled by an aluminum heat sink and not by a fan which is cause of unreliability, noise, and attracts dust,
- **Additional required assembly (extra paid):**
 - Telescopic trolley base,
 - Arm wall bracket,
 - Wall bracket,
 - Clamp for medical rails,
 - Universal clamp,
 - Simple clamp,
- Technical Passport,
- 3-year warranty

