



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **LID Bella LED light (New)**

- New,
- Made in France,
- A powerful lamp on a very sturdy, flexible arm, appreciated by professionals who want to perform minor procedures in their practice and who care about accurate rendering of skin and tissue color. The ergonomic handle allows precise and intuitive positioning,
- Dual light intensity option, 50% or 100%, allows modulation of light to suit all examinations,
- Specifications:
  - Light intensity: 1680 Lumens,
  - Light color: 4500 K (neutral white),
  - Color rendering index: IRC 95/100,
  - Electricity consumption: 13 W,
  - Class 1 medical device,
  - Lifetime of LEDs minimum 50,000 hours of operation,
  - Electrical protection: insulation class II,
  - Standards: IEC 60 601 -1/ -2-41 IEC 60 598-2-25,
  - Power cord with a length of 3 m,
  - Available colors: white,
- Dimensions :
  - Light head dimension: 160 x 130 mm,
  - Base dimension: 69 x 69 x 115 mm,
  - Length: flexible arm 700 mm,
  - Inter-axis spacing (4 holes in the lamp base): 54.5 mm,
- Weight: 2.9 kg,
- Light field diameter :
  - Distance 30 cm: 210,000 lux field illumination 15 cm
  - Distance 50 cm: 120,000 lux field illumination 20 cm,
  - Distance 100 cm: 33,000 lux field illumination 30 cm,
- Dual intensity switch optional at extra charge (on top of the lamp),
- Ergonomic handle allows precise placement of the light,
- Easy to clean and disinfect due to smooth surface and good quality materials,
- LEDs are implanted on a printed circuit cooled by an aluminum heat radiator, rather than by a fan, which is the cause of unreliability, noise and attracts dust,
- Additional mounting accessories are required for the lamp , your choice (at an extra charge):
  - Telescopic stand base 60-100 cm white metal, heavy duty
  - Telescopic stand base 60-100 cm white metal, heavy - version with brakes in 2 wheels
  - Telescopic stand base 90-130 cm white metal, heavy duty
  - Telescopic stand base 90-130 cm white metal, heavy - version with brakes in 2 wheels
  - Plain clamp - black
  - Clamp for medical rails with vertical lamp outlet
  - Universal clamp
  - Wall bracket
  - ARM wall bracket white (2x40cm)
  - ARM wall mount white (1x60cm)
  - Tabletop bracket ARM white (1x60cm)
  - Plate with metal tube in the middle  $\varnothing$ 16
  - Plate with metal tube on the edge  $\varnothing$ 16
  - Female socket with adaptor plate  $\varnothing$  16 mm x 27 mm
  - Adaptor plate male rod  $\varnothing$  15 mm x 30 mm,
- It has a valid Technical Passport issued,
- Warranty: 3 years,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

