



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## LED Lamp LID - Bella (New)

- New,
- Made in France,
- Watts: 17 W,
- Kelvins: 4000 K,
- Lumens: 1440 Lm,
- Light source: 12 LED,
- LED lifetime: 50 000 hours minimum,
- Color Rendering Index (IRC) 95,
- Power supply cable 3 m,
- Insulation Class II,
- Available colors : white
- Dimensions :
  - Head: 160 x 133 mm,
  - Base: 69 x 69 x 115 mm,
  - Weight: 2,3 kg,
  - Length: 700 mm,
- Diameter of illuminated field :
  - Distance 30 cm: 60000 lux, illuminated field: 5 cm
  - Distance 50 cm: 45000 lux, illuminated field: 10 cm,
  - Distance 100 cm: 15000 lux, illuminated field: 20 cm,
- Double intensity switch in option (on the top of the lamp),
- The ergonomic handle allows a precise and intuitive positioning,
- Easy cleaning and disinfection thanks to the smooth surfaces and quality materials,
- LED are implanted on a printed circuit cooled by an aluminum heat sink and not by a fan which is cause of unreliability, noise, and attracts dust,
- **Additional required assembly (extra paid):**
  - Telescopic trolley base,
  - Arm wall bracket,
  - Wall bracket,
  - Clamp for medical rails,
  - Universal clamp,
  - Simple clamp,
- Technical Passport,
- 3-year warranty

