



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Surgical Laser CO2 Surgery (New)

- New,
- Made in South Korea,
- Used in major specializations:
 - Dermatology,
 - Surgery,
 - Laryngology,
 - Proctology,
 - Gynaecology,
- Technical Specifications:
 - Wavelength: 10600 nm
 - Beam Mode Quality: TEM00,
 - Laser power output: 1 - 35 W,
 - Peak power: 315 W at 90 us,
 - Spot size: 120 - 300 μ m,
 - Repeat delay: 2 - 500 ms,
 - Aiming Beam: 3,5 mW, 650 nm,
 - Cooling System: liquid,
 - Ultra Puls Mode : 10, 12, 15- 50 mj,
 - Super Puls Mode: 55 - 100 mj,
 - Normal Puls Mode: 150 - 600 mj,
 - Continuous Wave Mode: 1 - 35 W,
 - Maximum power consumption: 1500 VA,
 - Power input: 240 V, 50/60 Hz,
 - Dimensions: 360 x 430 x 1200,
 - Weight: 45 kg,
- Technical Passport,
- 24 months warranty

